Precise Determination of the Mn-53/Mn-55 Isotopic Ratio in the ALH77250 Iron Meteorite

Shigekazu Yoneda¹, Takayuki Nagamine², Yasuji Oura², and Mitsuru Ebihara²

Department of Science and Engineering, National Science Museum,
 3-23-1 Hyakunin-cho, Shinjuku-ku, Tokyo 169-0073, Japan
 Department of Chemistry, Tokyo Metropolitan University,
 1-1 Minami-Osawa, Hachioji, Tokyo 192-0397, Japan

Abstract The natural radiogenic nuclide 53 Mn is found in iron meteorites. In order to determine the half life of 53 Mn, the 53 Mn/ 55 Mn isotopic ratios in the ALH77250, which has high 53 Mn concentration, were measured precisely by the thermal ionization mass spectrometer. The average of the isotopic ratios, when the measured ion currents of 55 Mn were greater than 1×10^{-12} A and the interference of 53 Cr to 53 Mn was less than 10%, was $(5.977\pm0.032)\times10^{-5}~(2\sigma_m)$. With combination of the X-ray measurements of 53 Mn, a new precise half life will be determined.

Key words: Mn-53, ALH77250, iron meteorite, isotopic ratio, half life

1. Introduction

Manganese-53 is one of the most useful natural radioactive nuclide for meteorite chronology. It is mainly produced in iron meteorites by cosmic rays with the reactions of ⁵⁶Fe(n, p3n), ⁵⁴Fe(n, 2n), or ⁵⁶Fe(p, a) (Honda et al.; 1961, Millard, 1965; Imamura et al., 1969), and called as a cosmogenic nuclide (e.g. Yabuki and Shima, 1989; Nishiizumi, 1987). But it is also known as one of extinct nuclides, which are once existed in the early solar system but already decayed due to short half lives compared with the age of the solar system $(4.56 \times 10^9 \text{ yr})$. The decay daughter isotope 53Cr is detected in the calcium, aluminum-rich inclusions in many carbonaceous chondrites and other meteorites, and it is used for formation age determination and decipher early history of the solar system (e.g. Birck and Allegre, 1988; Lugmair and Shukolyukov 1998). However, current value of the half life of 53Mn (3.7×10⁶ yr) is determined in 1971 (Honda and Imamura, 1971) and even the latest report is in 1972 (Wölfle et al., 1972). Thus we re-evaluate the half life of 53Mn and started to determine a new precise value. Here we report the measurements of the ⁵³Mn/⁵⁵Mn isotopic ratios which are of the order of 10⁻⁵ within 1% error, which is an important part to determine the half life.

2. Experiments

2.1 Sample Solutions

The Mn standard solution (Johnson-Matthey) was used for checking of the 55Mn ion beam at mass spectrometry. The method of the extraction of Mn from iron meteorites will be reported elsewhere and thus described here briefly. First, iron meteorites were dissolved in 1:2 HNO₃. Then, the solutions were dried up and converted to 7 M HCl solutions. With diisopropyl ether, Fe was extracted out. From the aqueous fraction, Mn was separated using anion and then cation exchange columns. The extracted Mn solutions were purified several times by the modified method of Honda and Imamura (1971). The method was tested by the Canyon Diablo iron meteorite and then applied to the ALH77250 iron meteorite, which has high concentration of 53Mn (493 dpm/kg; Setoguchi, 2001).




Fig. 1. A schematic diagram of a Micromass Sector 54–30 mass spectrometer. While ⁵⁵Mn is measured by the Faraday detector system, ⁵³Mn is measured by the Daly detector system through the electrostatic filter (ESF), which is operated in ion counting mode.

2.2 Instrumentation

A Micromass Sector 54-30 thermal ionization mass spectrometer at the National Science Museum was used for this study. A schematic diagram of the mass spectrometer is shown in Fig. 1. Seven Faraday cup collectors and a Daly detector system are placed at the end of the magnet and at the end of the electrostatic filter (ESF), respectively. While 55Mn is measured by the Faraday detector system, 53Mn is measured by the Daly detector system through the ESF. The Daly detector system is operated in the ion counting mode, which can count each 53Mn ion and can measure the ion beam current of as low as 1×10^{-18} A. Measuring 55Mn at a Faraday detector (up to 1× 10⁻¹⁰ A) and ⁵³Mn at the Daly detector, extreme isotope ratios of ⁵³Mn/⁵⁵Mn at the order of 10⁻⁵ can be measured within 1% analytical errors.

2.3 Mass Spectrometry of Mn

A rhenium ribbons were prebaked with 4.8 A filament current. Silica gel was placed on the ribbon in order to enhance the beam intensity of Mn and to reduce that of organic compounds. About $1 \mu g$ of Mn was loaded onto the ribbon and then

the sample was covered by the phosphoric acid. The ⁵³Mn/⁵⁵Mn ratios were determined via three sequences (Fig. 2) as in Yoneda *et al.* (1995). Sequence 2 measured ⁵³Mn at the Daly detector and ⁵⁵Mn at the H2 Faraday cup detector simultaneously. The sequences 1 and 3 were to determine the Cr interference on ⁵³Mn. The natural isotopic ratio of ⁵³Cr/⁵²Cr for the correction was assumed since contribution of spallation-produced ⁵³Cr is reported to be negligible (Shima and Honda, 1966).

Two separate measurements which used different filaments were performed. At each measurement, several sets of isotopic ratios in different filament currents were measured.

3. Results and Discussion

3.1 Faraday-Daly Gain Calibration

Just before the 53 Mn/ 55 Mn isotopic ratio measurements, the gain difference between Faraday and Daly collectors were measured by Faraday-Daly gain calibration. This calibration was performed by measuring the same 55 Mn ion beam of $\sim 10^{-13}$ A by Faraday and Daly collectors alter-




Fig. 2. The sequences for the ⁵³Mn/⁵⁵Mn measurements. At the sequence 2, ⁵³Mn and ⁵⁵Mn are measured simultaneously. The amount of the interference of ⁵³Cr on ⁵³Mn is calculated from the ⁵²Cr intensity at the sequence 1. The sequence 3 monitors Cr and other (such as Fe) interferences.

Table 1. Results of Faraday-Daly gain calibration. Measurement #1

Run	Filament current (A)	Mn-55 (×1E-13 A)	Daly/Faraday (×1E-14 A)	Standard error (%)
1	1.780	2.0	86.18	0.07
2	1.780	3.0	86.63	0.02
3	1.780	4.0	86.82	0.08
4	1.750	3.3	87.36	0.03
5	1.720	2.6	87.74	0.11
6	1.700	2.2	87.77	0.11
7	1.700	2.2	87.59	0.03
8	1.725	3.1	87.30	0.09
9	1.650	1.4	88.05	0.06

Mea	surement	#2
-----	----------	----

Run	Filament current (A)	Mn-55 (×1E-13 A)	Daly/Faraday (×1E-14 A)	Standard error (%)
1	1.810	2.2	84.09	0.04
2	1.810	2.6	84.36	0.04
3	1.810	2.9	84.63	0.01
4	1.814	3.4	85.15	0.02
5	1.814	3.5	85.37	0.03
6	1.800	3.1	85.95	0.04
7	1.790	2.8	85.97	0.01
8	1.780	2.5	86.17	0.03
9	1.770	2.3	86.24	0.02
10	1.790	2.8	85.98	0.02
11	1.800	3.2	85.78	0.03
12	1.810	3.4	85.71	0.03

Each run consists of 5 mesurement cycles.

nately. The results are shown in Table 1. The best values, 87.59 and 85.98, are used for the measurement #1 and #2, respectively. Using the ²³²Th beam, the Faraday-Daly gain factor variation was less than 0.5% during the 8-hour measurements (Yoneda *et al.*, 1995). The difference between the gain factors of two measurements is due to the slight difference of the settings of the mass spectrometer.

3.2 53Mn/55Mn Isotopic Ratio

The results of 2 separate measurements are listed in Table 2. The ion beam intensities of

⁵⁵Mn measured by the Faraday detector were from ~ 1 to 500×10^{-14} A, and those at mass 53 (mixture of ⁵³Mn and ⁵³Cr) measured by the Daly detector were from ~ 3 to 300×10^{-18} A, giving an isotopic ratio of $\sim 6\times 10^{-5}$. The correction of the ⁵³Cr interferences on ⁵³Mn calculated from the intensity of ⁵²Cr were normally small (less than 10% of the total intensity at the best condition). The average of each run, which is calculated from 10 ratios (except measurement #1 run 4) at the same filament current, are plotted against the filament currents in Fig. 3. The mass fractionation effect with the filament current, hence the

Table 2. Results of the Isotopic measurements of Mn-53/Mn-55. Measurement #1

	Filament current (A)	Cr-52	Mn-53+Cr-53 (×1E-14 A)	Mn-55 (×1E-14 A)	Cr-54+Fe-54 (×1E-14 A)	53Mn/55Mn (×1E-6)	1 Cr-53 (Mn-53+Cr-53) (%)	Avarage (×1E-6)	Standard error (%)
1	2.145	0.08887	0.01808	140.495	0.00359	56.959	55.74	58.922	0.914
		0.07236	0.01598	137.235	0.00327	56.654	51.35		
		0.05833	0.01434	134.240	0.00281	57.550	46.12		
		0.04553	0.01278	131.531	0.00223	57.910	40.40		
		0.03686	0.01210	128.562	0.00207	61.576	34.55		
		0.03035	0.01085	125.431	0.00188	59.078	31.72		
		0.02535	0.01019	122.470	0.00188	59.760	28.20		
		0.02213	0.00961	119.502	0.00158	59.427	26.11		
		0.01886	0.00933	117.100	0.00154	61.410	22.92		
		0.01657	0.00857	113.633	0.00144	58.901	21.92		
2	2.175	0.00496	0.00724	112.689	0.00103	59.216	7.78	61.029	0.673
		0.00447	0.00752	110.751	0.00101	63.318	6.75		
		0.00390	0.00709	106.987	0.00097	62.181	6.23		
		0.00359	0.00659	103.208	0.00090	59.917	6.18		
		0.00295	0.00655	99.399	0.00094	62.574	5.11		
		0.00261	0.00617	96.501	0.00089	60.861	4.80		
		0.00209	0.00597	94.206	0.00085	60.816	3.97		
		0.00187	0.00571	90.015	0.00068	61.045	3.72		
		0.00171	0.00535	86.148	0.00073	59.887	3.62		
		0.00151	0.00517	82.667	0.00060	60.470	3.31		
3	2.250	0.00157	0.00572	89.648	0.00061	61.832	3.12	59.420	1.098
		0.00145	0.00519	85.286	0.00058	58.955	3.17		
		0.00144	0.00502	81.283	0.00055	59.727	3.25		
		0.00133	0.00488	75.942	0.00053	62.318	3.08		
		0.00137	0.00465	72.769	0.00051	61.836	3.33		
		0.00137	0.00371	61.978	0.00039	57.391	4.18		
		0.00129	0.00429	70.011	0.00048	59.170	3.42		
		0.00133	0.00389	65.127	0.00041	57.441	3.89		
		0.00123	0.00394	64.153	0.00036	59.237	3.55		
		0.00120	0.00375	64.177	0.00040	56.288	3.63		
4	2.500	0.35233	0.04227	38.129	0.01008	60.751	94.52	39.714	19.957
		0.33806	0.03944	31.155	0.00929	35.566	97.19		
		0.30979	0.03573	26.663	0.00835	22.569	98.32		
		0.27936	0.03249	20.273	0.00779	39.968	97.51		

Table 2. Results of the Isotopic measurements of Mn-53/Mn-55 (continued). Measurement #2

Run	Filament current (A)	Cr-52	Mn-53+Cr-53 (×1E-14 A)	Mn-55 (×1E-14 A)	Cr-54+Fe-54 (×1E-14 A)	53Mn/55M1 (×1E-6)	Cr-53 (Mn-53+Cr-53) (%)	Avarage (×1E-6)	Standard error (%)
1	2.150	0.25660 0.24350 0.24020 0.22823 0.21353 0.20266 0.19171 0.18503 0.17610 0.16654	0.04312 0.04223 0.04063 0.03891 0.03820 0.03643 0.03438 0.03428 0.03383 0.03156	212.414 212.067 211.592 210.720 209.715 208.780 207.637 207.104 206.704 206.013	0.00808 0.00842 0.00871 0.00787 0.00726 0.00746 0.00692 0.00682 0.00640 0.00629	66.042 68.935 63.310 61.835 66.690 64.420 60.879 64.210 67.058 61.533	67.47 65.38 67.03 66.51 63.39 63.08 63.23 61.21 59.03 59.83	64.491	1.307
2	2.175	0.09483 0.08992 0.08805 0.08193 0.08004 0.08089 0.07681 0.07188 0.06563 0.06114	0.02674 0.02627 0.02592 0.02584 0.02501 0.02438 0.02381 0.02380 0.02262 0.02227	257.018 256.308 255.406 255.372 253.529 252.264 250.244 248.755 248.749 247.975	0.00000 0.00000 0.00000 0.00000 0.00000 0.00000 0.00000 0.00001	62.222 62.718 62.413 64.802 62.861 60.273 60.349 62.899 61.005 61.836	40.20 38.81 38.51 35.95 36.29 37.63 36.58 34.25 32.90 31.14	62.138	0.693
3	2.200	0.03909 0.03639 0.03519 0.03368 0.02956 0.03133 0.02873 0.02693 0.02571 0.02455	0.02313 0.02273 0.02278 0.02180 0.02208 0.02195 0.02180 0.02131 0.02098 0.02025	306.959 304.543 302.757 300.789 298.642 297.933 297.059 295.569 293.728 292.548	0.00000 0.00000 0.00000 0.00000 0.00000 0.00001 0.00000 0.00001 0.00000	60.910 61.099 62.073 59.788 62.703 61.737 62.406 61.783 61.490 59.719	19.17 18.15 17.51 17.52 15.18 16.19 14.95 14.33 13.90 13.74	61.371	0.520
4	2.250	0.02114 0.02152 0.01877 0.01652 0.01634 0.01511 0.01542 0.01453 0.01376 0.01287	0.02747 0.02762 0.02656 0.02634 0.02535 0.02482 0.02459 0.02486 0.02419 0.02384	424.372 419.943 415.254 408.418 390.641 384.896 382.489 385.081 379.212 371.105	0.00000 0.00000 0.00000 0.00000 0.00000 0.00000 0.00000 0.00001 0.00000 0.00000	59.083 59.972 58.844 59.916 60.146 60.037 59.722 60.283 59.670 60.311	8.73 8.83 8.01 7.11 7.31 6.90 7.11 6.63 6.45 6.12	59.798	0.259
5	2.300	0.04754 0.04209 0.03661 0.03743 0.03189 0.03064 0.02937 0.02573 0.02433 0.02343	0.03746 0.03624 0.03404 0.03394 0.03144 0.03068 0.02979 0.02939 0.02855 0.02806	520.598 512.763 494.381 479.104 460.116 454.062 439.623 434.910 424.409 418.994	0.00000 0.00000 0.00000 0.00000 0.00000 0.00000 0.00001 0.00001 0.00000	61.610 61.364 60.459 61.972 60.468 59.912 60.179 60.857 60.772 60.637	14.39 13.17 12.20 12.51 11.50 11.33 11.18 9.93 9.66 9.47	60.823	0.337

Table 2. Results of the Isotopic measurements of Mn-53/Mn-55 (continued). Measurement #2 (continued)

6 2.350 0.01850 0.02883 440.098 0.00000 60.752 0.01544 0.02637 414.709 0.00000 59.364 0.01374 0.02547 397.589 0.00001 60.139 0.01267 0.02440 386.342 0.00000 59.446 0.01138 0.02286 371.609 0.00000 58.044 0.01009 0.02239 357.878 0.00001 59.359 0.00907 0.02150 346.232 0.00001 59.137 0.00787 0.02099 338.677 0.00001 59.335 0.00752 0.01993 329.389 0.00000 57.917 0.00720 0.01934 314.951 -0.00001 58.812	7.28 6.64 6.12 5.89 5.65 5.11 4.78 4.25 4.28 4.22 2.74	59.230	0.455
0.01544 0.02637 414.709 0.00000 59.364 0.01374 0.02547 397.589 0.00001 60.139 0.01267 0.02440 386.342 0.00000 59.446 0.01138 0.02286 371.609 0.00000 58.044 0.01009 0.02239 357.878 0.00001 59.359 0.00907 0.02150 346.232 0.00001 59.137 0.00787 0.02099 338.677 0.00001 59.335 0.00752 0.01993 329.389 0.00000 57.917	6.12 5.89 5.65 5.11 4.78 4.25 4.28 4.22		
0.01374 0.02547 397.589 0.00001 60.139 0.01267 0.02440 386.342 0.00000 59.446 0.01138 0.02286 371.609 0.00000 58.044 0.01009 0.02239 357.878 0.00001 59.359 0.00907 0.02150 346.232 0.00001 59.137 0.00787 0.02099 338.677 0.00001 59.335 0.00752 0.01993 329.389 0.00000 57.917	5.89 5.65 5.11 4.78 4.25 4.28 4.22		
0.01138 0.02286 371.609 0.00000 58.044 0.01009 0.02239 357.878 0.00001 59.359 0.00907 0.02150 346.232 0.00001 59.137 0.00787 0.02099 338.677 0.00001 59.335 0.00752 0.01993 329.389 0.00000 57.917	5.65 5.11 4.78 4.25 4.28 4.22		
0.01009 0.02239 357.878 0.00001 59.359 0.00907 0.02150 346.232 0.00001 59.137 0.00787 0.02099 338.677 0.00001 59.335 0.00752 0.01993 329.389 0.00000 57.917	5.11 4.78 4.25 4.28 4.22		
0.00907 0.02150 346.232 0.00001 59.137 0.00787 0.02099 338.677 0.00001 59.335 0.00752 0.01993 329.389 0.00000 57.917	4.78 4.25 4.28 4.22		
0.00787 0.02099 338.677 0.00001 59.335 0.00752 0.01993 329.389 0.00000 57.917	4.25 4.28 4.22		
0.00752 0.01993 329.389 0.00000 57.917	4.28 4.22		
***************************************	4.22		
0.00720 0.01024 214.051 -0.00001 59.912			
0.00720 0.01934 314.951 -0.00001 58.812	2.74		
7 2.400 0.00461 0.01912 315.103 0.00000 59.015		59.482	0.495
0.00398	2.44		
0.00351 0.01753 285.610 0.00001 59.974	2.27		
0.00316 0.01684 275.183 0.00001 59.894	2.13		
0.00307	2.24		
0.00264 0.01350 225.649 0.00000 58.511	2.22		
0.00236	2.15		
0.00217	2.09		
0.00195 0.01067 174.387 0.00001 59.905	2.07		
0.00167 0.00966 162.268 -0.00001 58.372	1.96		
8 2.500 0.00292 0.00562 84.853 0.00052 62.353	5.89	59.847	0.774
0.00290 0.00527 81.597 0.00059 60.551	6.23		
0.00272	6.39		
0.00245	6.39		
0.00247	6.43		
0.00203 0.00398 62.688 0.00036 59.762	5.78		
0.00221 0.00398 61.127 0.00034 61.047	6.29		
0.00204 0.00356 57.812 0.00036 57.590	6.49		
0.00195	7.78 9.07		
0.00185 0.00232 34.918 0.00023 60.395		5 0.040	0.555
9 2.700 0.00636 0.01125 175.443 0.00073 60.024	6.41	59.848	0.755
0.00614	6.66		
0.00575 0.01017 165.204 0.00058 57.589	6.41		
0.00543	6.28		
0.00495 0.00915 143.485 0.00039 59.825	6.14		
0.00455 0.00819 126.719 0.00032 60.543	6.30		
0.00390	6.51		
0.00360 0.00650 100.591 0.00024 60.549	6.28		
0.00311 0.00544 87.460 0.00020 58.147 0.00294 0.00558 83.933 0.00016 62.500	6.49 5.97		
	21.98	57.399	4.866
10 2.850 0.00325 0.00167 20.646 0.00011 63.265 0.00345 0.00140 15.918 0.00012 63.341	27.93	31.377	1.000
	38.81		
0.00316 0.00092 11.764 0.00013 48.016 0.00303 0.00085 8.757 0.00010 57.337	40.63		
0.00305 0.00085 8.757 0.00010 57.337 0.00305 0.00075 6.247 0.00010 63.893	46.45		
0.00305 0.00075 0.247 0.00010 03.093	52.49		
0.00284 0.00050 3.011 0.00009 60.684	63.78		
	74.75		
0.00287 0.00044 2.008 0.00011 54.792 0.00271 0.00038 1.363 0.00008 52.481	81.09		
0.00271 0.00038 1.303 0.00008 32.401 0.00251 0.00032 0.898 0.00009 40.099	88.79		
0.00251 0.00032 0.898 0.00009 40.099	00.77		




Fig. 3. The results of the 53 Mn/ 55 Mn ratios of measurement #1 (open circles) and measurement #2 (solid squares) plotted against the filament currents. Error bars represent the standard errors ($2\sigma_{\rm m}$).

temperature, are not seen during the measurements. For the calculation of the total average of the isotopic ratios, only the data, when the measured ion currents of ^{55}Mn were greater than $1\times 10^{-12}\,\text{A}$ and the interference of ^{53}Cr to ^{53}Mn was less than 10%, were used. The obtained average was $(5.977\pm0.032)\times10^{-5}$ (errors are $2\,\sigma_{\text{m}}$).

3.3 Half life of 53Mn

The half life of 53 Mn $(T_{1/2})$ can be calculated from the equation (1):

$$T_{1/2} = \frac{\ln 2 \times N(^{53} \text{Mn})}{A(^{53} \text{Mn})} \tag{1}$$

where, $N(^{53}\text{Mn})$ is number of ^{53}Mn atoms and $A(^{53}\text{Mn})$ is decay rate. $N(^{53}\text{Mn})$ is can be calculated from the isotopic ratio $^{53}\text{Mn}/^{55}\text{Mn}$ and the content of Mn in the sample. $A(^{53}\text{Mn})$ is determined by X-ray countings with a LEPS counter. Preliminary results of these measurements were reported by Oura *et al.* (2002). With 0.09289 Bq and 19.14 μ g Mn in the sample, a new half life of

 53 Mn, 3.0×10^6 y, was calculated. This new value is about 20% shorter than the current value 3.7×10^6 y and the relative analytical errors of the new value, 5%, is smaller than those of the current value, 10%. However, the errors are still large comparing with the isotopic ratio measurements for meteorite chronology. Additional X-ray measurements and efforts to reduce uncertainties are in progress. Details will be published elsewhere.

Acknowledgments

We are thankful to Prof. M. Honda of Nihon University and Prof. M. Imamura of the Institute for Nuclear Studies, the University of Tokyo, for providing us the ⁵³Mn standard solution and helpful comments, and Prof. H. Kojima of the National Institute of Polar Research, for providing us the ALH77250 iron meteorite. This work was supported in part by a Grant-in-Aid for Scientific Research of Ministry of Education, Science and Culture, Japan (No. 05453012 to M.E.).

References

- Birck and Allègre, 1988. Manganese-chromium isotopic systematics and development of the early solar system. *Nature*, **331**, 579–584.
- Honda M. and M. Imamura, 1971. Half-life of Mn⁵³. Physical Review, 4: 1182–1188.
- Honda, M., J. P. Shedlovsky and J. R. Arnold, 1961. Radioactive species produced by cosmic rays in iron meteorites. *Geochim. Cosmochim. Acta*, 22: 133–154.
- Imamura, M., H. Matsuda, K. Horie and M. Honda, 1969.
 Application of neutron activation method for ⁵³Mn in meteoritic iron. *Earth Planet. Sci. Lett.*, 6: 165–172.
- Lugmair G. W. and A. Shukolyukov 1998. Early solar system timescales according to ⁵³Mn–⁵³Cr systematics. *Geochim. Cosmochim. Acta*, **62**, 2863–2886.
- Millard, H. T. Jr., 1965. Thermal neutron activation: Measurement of cross section for manganese-53. *Science*, 147: 503–504.
- Nishiizumi, K., 1987. ⁵³Mn, ²⁶Al, ¹⁰Be and ³⁶Cl in meteorites: Data compilation. *Nucl. Tracks Radiat. Meas.*, **13**: 209–273.

- Oura, Y., T. Nagamine, M. Ebihara, S. Yoneda and M. Honda, 2002. Measurement of half life of ⁵³Mn (III) (abstract), Annual Meeting of the Radiochemical Society of Japan 2002.
- Setoguchi, M., 2001. Master theses, Tokyo Metropolitan University.
- Shima, M. and M. Honda, 1966. Distribution of spallation produced chromium between alloys in iron meteorites. *Earth Planet. Sci. Lett.*, 1: 65–74.
- Wölfle, R., W. Herr and U. Herpers 1972. Determination of the activation cross section σ_{therm} and of the halflife of ⁵³Mn by a 587 days reactor irradiation. *Radiochim. Acta*, **18**, 207–211.
- Yabuki, S. and M. Shima, 1989. Cosmogenic ²⁶Al and ⁵³Mn in Japanese chondrites. *Bull. Natn. Sci. Mus.*, Ser. F. 12: 1–14.
- Yoneda, S., M. Shima, H. Hidaka and M. Ebihara, 1995. An attempt of isotopic abundance determination of natural radioactive nuclides, ⁵³Mn and ²³⁰Th, by thermal ionization mass spectrometry. *Bull. Natn. Sci. Mus., Ser. E*, **18**: 1–8.